



HOW TO RECONSTITUTE BAC WATER

Step-by-Step Guide

1 Gather Supplies

- BAC water vial,
- sterile powder vial,
- syringe,
- alcohol prep pads



2 Sanitize Vials

- Use alcohol prep pad to clean the tops of both vials



3 Inject BAC Water

- Inject all 10 mL of BAC water into the powder vial.
- -Do not shake



4 Mix Gently

- Gently swirl the vial until the powder fully dissolves.



5 Store Refrigerated

- Once reconstituted, keep the vial refrigerated



Once peptide is reconstituted, it must be refrigerated.

For Research & Educational Use Only
Not for human or veterinary use